Symbols for Warnings, Cautionary Measures and Notes

The following conventions are used in this manual:

WARNING!

This symbol advises the user of serious danger for the patient and the user.

CAUTION!

The symbol informs the user that particular care is required for safe and efficient operation of the system.

IMPORTANT NOTE

This symbol provides the user with useful or additional information.

Product Description

The WorldPRP hourglass tube is a single-use, sterile device that is used for separating blood. The components are outlined in Figure 1.


Figure 1: Components of the WorldPRP tube

Intended Use

The WorldPRP tube is indicated for platelet rich plasma (PRP) separation from whole blood via a centrifuge.

Contraindications:

The WorldPRP is not to be used under the following contraindications:

1. Critical thrombocytopenia (low platelet count)
2. Hypofibrinogenemia
3. Haemodynamic instability (collapse)
4. Sepsis (infection)
5. Acute and chronic infections
6. Chronic liver disease
7. Anti-coagulation therapy (warfarin, dabigatran, heparin)

WARNING!

• ONLY inject the PRP mixture into the patient from whom the PRP was derived

- Process only one patient's blood per disposable WorldPRP device
- Do NOT use for direct intravenous infusion
- Do NOT re-use the used WorldPRP device
- Do NOT re-sterilize
- Do NOT use expired products
- Ensure patient is aware of general risks associated with treatment and possible adverse effects

CAUTION!

- ONLY use certified blood collecting equipment
- Always store the device in clean state
- Ensure device is not contaminated
- Ensure Identify tag is affixed to the WorldPRP (containing donor's name or registration No.) to avoid contamination
- Ensure drawn blood is not exposed to other people
- Discard disposable device in biohazard container after use
- PRP treatment is to be carried out by trained Health Professional

How to use:

Prior to Use

1. Check the expiry date of the product.
2. Check for any damage on the packaging and product.

WARNING!

- Do NOT use if product date is expired
- Do NOT use if inner packaging seal is broken as it is no longer sterile
- Do NOT use if damaged

CAUTION!

- Read all instructions prior to use

Instructions (based on a 23ml whole blood sample + 2ml anticoagulant)

1. Ensure the centrifuge is on a stable surface, power cord attached, and the unit is turned on.
2. With gloved hands, unpack the WorldPRP tube from the sterile packaging.
3. Ensure top chamber is open to bottom chamber (to open turn top collar counterclockwise one and a quarter turns).
4. Ensure bottom chamber is open to 15ml mark (to adjust, gently rotate bottom collar clockwise until the top of the internal piston is at the 15ml mark).
5. Open outer port by lifting cap. Transfer the whole blood sample into the WorldPRP tube using the syringe (without a needle) and inserting into outer port.
 - a. Do not push the plunger with force, lean the tube at less than a 45° angle and gently transfer the whole blood. Reseal the outer port after the transfer by closing cap.

WARNING!

- Do NOT push plunger with strong force when injecting into the WorldPRP tube

6. Using a digital scale, ensure that the weight of the WorldPRP tube and the counter weight are equal (\pm 5g). Place the fully prepared WorldPRP tube and counter weight tube into the centrifuge, aligned across from each other. Close the lid and lock the centrifuge.
7. Place the WorldPRP tube and counter balance into the centrifuge. Set the centrifuge to the following settings: Rotation speed 3,000 rpm, Operation time: 3 minutes.

*For additional centrifuge information please refer to the user manual.

IMPORTANT NOTE

- Follow the instructions provided with the centrifuge regarding how to set the required settings

8. Once the cycle has stopped, unlock the lid and slowly remove the hourglass tube.
9. The red blood cells (RBC) will be located in the bottom chamber. The buffy coat (BC) sits immediately on the top of the RBC in the bottom chamber of the device. Raise the RBC and BC by turning the bottom collar counterclockwise until the buffy coat has reached the neck of the hourglass.
10. After ensuring the buffy coat has been positioned into the top chamber turn the top collar clockwise, ensuring inner port has sealed red blood cells in the bottom chamber.
11. Collect PRP from outer port using a syringe (no needle).

a. Possible collection options

OPTION 1: Extra-concentrated

For a more concentrated version, immediately after sealing the RBCs in the bottom chamber tilt the hourglass towards you 45° with outer port at 6 o'clock position and aspirate half of the contents. The aspirated contents are the PRP & buffy coat. This will leave the platelet poor plasma in the hourglass.

OPTION 2: Higher volume

Gently swirl the WorldPRP tube to mix the platelets, plasma, and BC and then aspirate the PRP through outer port.

12. Discard all biological materials as per local regulations.

Manufacturer Information

LUVO^{MD} 125 Fleming Drive
Cambridge, ON
N1T 2B8

For any questions, contact LUVO Medical authorized representative.

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